

Filing System AVS-1050/750



Stable construction made of powder coated steel (RAL 7035, grey) with 11 drawers with stopper. Cover plate of the drawers in blue color (RAL 5012) with handle bar made of aluminium and covered label.

Different models

| | | |
|---------------|--|----------------------|
| Type: 750-OT | Filing system for 113,300 glass slides | 800 x 750 x 1200 mm |
| Type: 1050-OT | Filing system for 169,400 glass slides | 1050 x 750 x 1200 mm |
| Type: 750-P | Filing system for 30,800 paraffin blocks | 800 x 750 x 1200 mm |
| Type: 1050-P | Filing system for 46,200 paraffin blocks | 1050 x 750 x 1200 mm |

Plastic trays for glass slides and paraffin blocks

| | |
|------------------|------------------|
| Dimension W/D/H: | 64 x 310 x 40 mm |
| Color: | RAL 7011 |
| Article-No.: | 090.500.100 |

Technical details see back page

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de



Filing System AVS-1050/750

Technical Data:

Filing system for glass slides – AVS-1050/750-OT

| | |
|--------------------|------------------------------|
| Blocks per drawer: | approx. 10,300 glass slides |
| Total amount: | approx. 113,300 glass slides |
| Dead weight: | approx. 205 kg |
| Overall weight: | approx.. 1,200 kg |
| Dimensions W/D/H: | 800 x 750 x 1200 mm |
| Type: | AVS-750-10-OT |
| Article-No.: | 090.500.010 |

| | |
|--------------------|------------------------------|
| Blocks per drawer: | approx. 15,400 glass slides |
| Total amount: | approx. 169,400 glass slides |
| Dead weight: | approx. 240 kg |
| Overall weight: | approx.. 1,750 kg |
| Dimensions W/D/H: | 1050 x 750 x 1200 mm |
| Type: | AVS-1050-10-OT |
| Article-No.: | 090.500.020 |

Filing system for paraffin blocks – AVS-1050/750-P

| | |
|--------------------|--------------------------------|
| Blocks per drawer: | approx. 2,800 paraffin blocks |
| Total amount: | approx. 30,800 paraffin blocks |
| Dead weight: | approx. 205kg |
| Overall weight: | approx. 590 kg |
| Dimensions W/D/H: | 800 x 750 x 1200 mm |
| Type: | AVS-750-10-P |
| Article-No.: | 090.500.030 |

| | |
|--------------------|--------------------------------|
| Blocks per drawer: | approx. 4,200 paraffin blocks |
| Total amount: | approx. 46,200 paraffin blocks |
| Dead weight: | approx. 240kg |
| Overall weight: | approx.. 1,050 kg |
| Dimensions W/D/H: | 1050 x 750 x 1200 mm |
| Type: | AVS-1050-10-P |
| Article-No.: | 090.500.040 |

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de

