


## Trolley Made of Stainless Steel RW-4 – RW-15 for Preparation Cabinets



The trolley is completely made of a stainless steel and is designed for storing tissue test glass containers. Consisting of 4 wheels  $d=125$  mm, 2 of them lockable and an anti-collision buffer.

Stainless steel containers are to be ordered separately and are available in: 20 mm, 40 mm and 65 mm depth

### **Technical Data:**

Dimensions W/D/H:	715 x 535 x 1039 mm
Numbers of shelves:	4
Type:	RW-4
Article No.:	090.020.225

*Additional models see back page*

**KUGEL medical GmbH & Co. KG**  
Hermann-Köhl-Str. 2A  
DE - 93049 Regensburg / Germany

Phone +49/941/208648-0  
Fax +49/941/208648-29  
[www.KUGEL-medical.de](http://www.KUGEL-medical.de)



**KUGEL**  
medical  
■ ■ ■ ■

# Trolley Made of Stainless Steel RW-4 – RW-15 for Preparation Cabinets

## **Model RW-6:**

Dimensions W/D/H: 715 x 535 x 1670 mm  
Numbers of shelves: 6  
Type: RW-6  
Article No.: 090.020.220

## **Model RW-9:**

Dimensions W/D/H: 715 x 535 x 1785 mm  
Numbers of shelves: 9  
Type: RW-9  
Article No.: 090.020.229

## **Model RW-15:**

Dimensions W/D/H: 715 x 535 x 1800 mm  
Numbers of shelves: 15  
Type: RW-15  
Article No.: 090.020.227

**KUGEL medical GmbH & Co. KG**  
Hermann-Köhl-Str. 2A  
DE – 93049 Regensburg / Germany

Phone +49/941/208648-0  
Fax +49/941/208648-29  
[www.KUGEL-medical.de](http://www.KUGEL-medical.de)



**KUGEL**  
*medical*  
■ ■ ■ ■