

Autopsy Table ST 10/100

Vertically Adjustable and Rotatable



Model ST 10/150 with additional tap beside the basin

The autopsy table is mounted on a center support column and can be rotated to the right and the left by approx. 180°; moreover, it is height-adjustable from 750 to 1000 mm. The autopsy table is completely made of stainless steel (material no. 1.4301- ANSI 304). The work surface features a high, profiled surround and a seamlessly welded basin, 400 x 500 x 200 mm, with a tube overflow and a sieve insert is installed at one end of the table. Additionally, the table top is equipped with a waterproof sound-proofing and features a negative inclination arranged diagonally to the drain. All power connections are built into the base and can be easily accessed by means of an inspection opening.

The Autopsy Table ST 10/100 model includes the following features:

- High-speed clamping grip to fix the rotating movement
- Knee-operated mixing tap for cold and warm water at the center column
- 3 m shower hose with hand sprinkler
- 2 splash proof electrical outlets
- UP / DOWN button for height adjustment by means of an integrated electro mechanical drive
- Electrical and water connections are preinstalled
- The safety requirements for drinking water as per DIN 1988 part 4 and EN DIN 1717 are to be provided by the customer

Technical Data:

| | |
|--------------------|----------------------------|
| Dimensions W/D/H: | 2600 x 850 x 750 / 1000 mm |
| Voltage: | 230 Volt / 1 Ph / 50 Hz |
| Carrying capacity: | 300 kg |
| Type: | ST 10/100 |
| Article No.: | 010.010.100 |

Alternatively available with a carrying capacity of 450 kg at an extra charge – Article No.: 010.010.858

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de

