

Electronic Organ Scale with Weighing Bowl EOW-30



The completely electronic scale for precise weighing has a large LCD display, height 20 mm. The weight is switchable from kg to lbs. The scale has an automatic zero setting, tare function and hold function. The stainless steel bowl can be removed for cleaning. Delivery includes batteries and main supply.

Technical Data:

Dimensions W/D/H:	535 x 241 x 330 mm
Weighing range:	0 g – 30 kg
Accuracy:	1 gram
Type:	EOW-30
Article No.:	090.010.095

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de



KUGEL
medical
■ ■ ■ ■